

Newsletter 2022

1 Temperature Measurements

2 Tube Rotator

3 Tips & Tricks

→ Aspiration Pump Syringe Purge

4 Publication Spotlight

→ Development of a Potent High-Affinity Human Therapeutic Antibody

5 Ask the Inventor

→ Linear Range

6 Geek Corner

→ Autosampler Calibration Wizard



1

Temperature Measurements

With the release of the *TC1000*, KinExA users are now able to run experiments at temperatures ranging from 4°C to 40°C. While useful for characterizing biological interactions, there are some things to consider when changing temperature.

Equilibrium & K_d shifts

Samples need to equilibrate at the temperature at which they are being evaluated. This means that samples prepared on the bench need to reach the temperature inside the *TC1000* and then reach binding equilibrium at that temperature. Small samples reach temperature more quickly than large samples. Make sure to give samples enough time to reach temperature in addition to equilibration time when calculating incubation time needed.

K_d values are temperature dependent. We have routinely seen K_d 's get weaker with higher temperatures but that is not always the case. If you have a system that you measured at room temperature, do not expect the same results at other temperatures. If the K_d becomes tighter, constant binding partner (CBP) and titrant concentrations may need to shift down to accurately define the K_d . Conversely, if the K_d becomes weaker, you may need to use higher concentrations. Note that CBP concentrations need to stay in linear range and if working with weaker K_d 's, make sure you are still in the KinExA Mode. Refer to *TN224 Linear Range* and *TN221 KinExA Mode* for more information.

Condensation

Switching from 4°C to warmer temperatures can cause condensation in the instrument. Normally, if the *TC1000* remains closed until the new temperature is reached, condensation will not be an issue. Note: the buffer in the reservoir will take much longer to reach the new temperature. To ensure condensation does not affect the instrument, and that buffer is at the right temperature, change the temperature at the end of the day. The instrument and buffer will be ready the next morning.

Evaporation

Small Volume Tubes (P/N: [111642](#)) are often used for small sample volumes and to reduce dead volume needed. According to *TN206 Minimum Sample Volume*, as little as 10 μ L is sufficient for dead volume for an experiment duration of 5 hours at room temperature. That will not be the case at 37°C; evaporation happens more quickly so additional volume should be added. While we recommend a minimum of 50 μ L dead volume, less can be evaluated and potentially used.

Samples should be covered during equilibration. This is true at any temperature. *Foil Tape* (P/N: [JVCCAF20](#)) or scotch tape are good options, but do not use Parafilm. When the experiment starts, the sample tip will pierce the tape to access the sample. As the experiment progresses, there is a hole that allows sample to evaporate. If that is a concern due to small samples, consider preparing the titration series and then splitting them in half to aliquot into additional tubes. Cover all samples and adjust the sample set to sample each tube once.

Temperature Tolerant Parts

The *XP* and *XCalibur* Pumps used on 3X00 and 4000 KinExA models come standard with plastic *XP* or *XCalibur* Valves. In order to avoid the temperature effects of plastic pieces, *Ceramic XP Valves* (P/N: [310136](#)) or *Ceramic Xcalibur Valves* (P/N: [310131](#)) will need to be used. Syringe barrels that contain a plastic piston should also be replaced with *Temperature Tolerant Syringe Barrels* (P/N [310510](#)). It is important that these *Valves* are in good working condition to prevent leaks that may lead to bubble and sample draw issues.

Note: These parts are included with the purchase of a *TC1000*.



TC1000 with KinExA 4000

2

Tube Rotator

Our tube rotator was designed to optimize KinExA experiment preparation. It employs gentle and effective mixing with an adjustable angle that optimizes mixing for any sample. The compact design is small enough to fit on a crowded bench, in an incubator, or under the hood. Purchase includes the base unit, power cord, one tube roller mixing rotor, and one end-over-end mixing rotor. Rotors can hold up to ten 50 mL tubes, ten 5 or 15 mL tubes, and up to 30 micro volume tubes.

Rotating Speed	10 RPM @ 60Hz 8 RPM @ 50 Hz
Wattage	4W
Voltage	12V
Mixing Angle	0 to 45
Height	11 in 279 mm
Length	9.5 in 241 mm
Width	5.5 in 139 mm



3

Tips & Tricks

Aspiration Pump Syringe Purge

Sapidyne's latest software release, KinExA Pro version 4.5, includes a new feature called the *Aspiration Pump Syringe Purge*. As the name suggests, this feature allows you to "purge" your instrument's syringe barrel of those little bubbles that sometimes cling to the top of the syringe. Keeping air out of the system helps in getting clean and accurate data.

To purge the syringe open the Instrument tab, click on the "Syringe Purge" button, then click "Start". If you are using a 3200 model, the software will ask you to loosen the Luer cap on the degasser module until the vacuum pump engages, then retighten the cap (see **Figure 1**). The KinExA 4000 can engage the vacuum pump through software alone, but the 3200 needs you to release some pressure to engage the vacuum pump. Afterwards, both models will automatically perform the syringe purge and inform you when it is complete (about 2 minutes). This feature only exists on KinExA Pro software version 4.5 and newer.




Figure 1. Degasser module on KinExA 3200 model with Luer Cap disconnected.

Publication Spotlight



Development of a potent high-affinity human therapeutic antibody via novel application of recombination signal sequence-based affinity maturation

Agnieszka Kielczewska, Igor D'Angelo, Maria Sheena Amador, Tina Wang, Athena Sudom, Xiaoshan Min, Palaniswami Rathanaswami, Craig Pigott, Ian N. Foltz

DOI: <https://doi.org/10.1016/j.jbc.2021.101533>

In a recent publication from Amgen, researchers used mammalian recombination signal sequence (RSS) as an affinity maturation method to improve the K_d of an antibody-antigen interaction. This method was able to produce lead antibodies ranging from 12 pM to 34 fM (**Figure 2**). Measuring interactions this tight can be challenging, even with the KinExA instrument.

The scientists at Amgen were able to achieve this by using multi curve analysis, and paying attention to the ratio (CBP concentration to K_d) of the binding curves. As discussed in a previous publication¹, it is important that the active binding site concentration be near or below the actual K_d for accurate K_d measurements. A CBP concentration as low as 50 fM was able to be measured, allowing for defined 95% confidence intervals. For the data shown in the **Fig. 2**, a 1:1 binding model was used to fit theory curves to the data. Each curve shows a good fit both in terms of shape and spacing of the curves. The 34 fM K_d measured here is the tightest antibody K_d published to date using the KinExA Instrument.




Figure 2. Affinity Determination of antibody variants: KinExA method affinity determination of lead antibodies to IL-13.

1. P. Rathanaswami et al. / Biochemical and Biophysical Research Communications 334 (2005) 1004–1013



5

Ask the Inventor

Why do measurements need to be performed with signals in the linear range? How will equilibrium & kinetic data be affected outside of the linear range?



When the KinExA Pro analysis converts binding signals (volts) into free Constant Binding Partner (CBP) concentrations (molar), it does so in part by assuming the signal is directly proportional to the free CBP concentration. Signals are generated by binding free CBP to Titrant coated solid phase and the linearity of response is based on an assumption that the binding is first order, meaning the fraction of solid phase Titrant that is bound is an insignificant fraction of the total solid phase Titrant. The KinExA linearity test (see *TechNote 224 Linear Range* for details) is a method validating the assumption, and provides guidance for the maximum CBP concentration that should be used with the specific beads and timing tested. As part of the Linear Range test, a hyperbolic equation (**Equation 1**) is fitted to measured data to get estimates for the A and B parameters.

$$\text{Measured Signal} = \left(\frac{A \times [\text{Free CBP}]}{[B] + [\text{Free CBP}]} \right) + \text{NSB}$$

Equation 1.

A = maximum saturation signal
 B = concentration that produces half of the max signal
 NSB = non specific binding signal

The maximum CBP concentration Sapidyne recommends is $0.2 \times B$ which is illustrated in **Figure 3**.

My guess is that none of my readers will object that Sapidyne's recommendation is not linear enough but that some will think it's too conservative and higher CBP concentrations could be used without problems. I acknowledge that $0.2 \times B$ is a judgment call on our part and not a hard and fast rule. With that in mind, I constructed a series of simulations to give you a sense of how much error you can expect in your K_d and k_{on} determinations for a variety of B factors ($=CBP/B$) should you choose to operate outside our guidance. I simulated a dual curve system with a K_d of 1 nM and CBP concentrations of 5 and 50 nM. I looked at B factors ranging from 0.01 to 5 which I implemented by changing B rather than changing the CBP. For each B factor I simulated hyperbolic signals, then added pseudorandom noise (3% of Sig100) and then solved for K_d using the standard KinExA Pro analysis. I repeated with a minimum of 50 and up to 1000 different random noise sets and then averaged the K_d value and the confidence intervals. The results are shown in **Figure 4** (equilibrium dual curve) and **Figure 5** (kinetics direct).

Note

The effect of the hyperbolic response is that signals are higher than they would be if the response was linear. Curves are shifted toward the right and made steeper the larger the CBP concentration is relative to B (See *Tech Note 224 Linear Range* for some examples). This means that curves measured with ratios that would normally easily resolve the K_d , for example a ratio of 20, may be analyzed as much higher ratios and fail to resolve the K_d . I saw this in simulations using a 200 pM low curve either alone or with a 20 nM high curve.




Figure 3. Typical hyperbolic response curve.





Figure 4. The factor of error in reported K_d value vs. the CBP concentration expressed as a fraction of B, the concentration of half maximum for the hyperbolic response. As an example, if the CBP concentration is equal to B (1 on the x axis) you can expect a factor of error between 1 (no error) and approximately 3x. The reported K_d 's are generally weaker than actual. The yellow circle represents Sapidyne's recommended range of operation, $CBP/B \leq 0.2$.

Figure 5. The factor of error in k_{on} as a function of the CBP expressed as a fraction of B. For example, if the CBP is equal to $2*B$ (2 on the x axis) the expected range for the reported k_{on} is between approximately 0.4 and 0.65 times the actual k_{on} . Reported k_{on} values are expected to be slower than actual. The yellow circle represents Sapidyne's recommended range of operation, $CBP/B \leq 0.2$.



Conclusions

Sapidyne's guidance is to operate at a concentration of $0.2*B$ or less. As you can see in **Figures 4 and 5** it is a case of the "lower the better". In many cases this means you are operating with Sig100 levels of less than a volt (itself a rule of thumb at Sapidyne). Many users prefer large signals but studies conducted at Sapidyne (see *TechNote 228 "Signal vs Signal/Noise"*) show that raising signals above 0.5 volts is not productive in terms of K_d accuracy.

Finally, if an experiment must be conducted using concentrations with a higher percentage of the B factor the error can be corrected by submitting the sanitized .kxp file along with the hyperbolic signal test to a Sapidyne representative.



KinExA Training Videos, How To Guides, & More
www.sapidyne.com/resources

6

Geek Corner

Autosampler Calibration Wizard



To achieve proper sampling during experiments and to avoid collisions with the sipper tip, the autosampler must be calibrated. In software versions 4.5.X and higher, improvements have been made that make this process easier. Previously, users were left to define calibrations independently and through numerical grid positioning. To calibrate the autosampler, navigate to the tools bar, drop to calibration and designate Autosampler. A screen, shown in **Figure 6**, will appear with lists of everything that can be calibrated. Positions that have saved calibrations will be marked by a green check and un-calibrated positions will have an empty circle (**Figure 7**).



Figure 7. Check mark indicates calibration status.

To calibrate a specific position, click on it and the calibration wizard (**Figure 8**) will guide you through the calibration process. Follow the directions on the calibration wizard pop-up, and click next when completed.

What makes the new calibrations significantly easier than previous software versions, is the ability to direct the sipper tip using arrows and also to change the increment in which it moves. Using the green arrows,




Figure 6. Autosampler calibration page.




Figure 8. Calibration wizard.

the sipper tip can be directed towards the center of the position of interest. The bottom most arrow will direct the sipper tip towards the user with all others following this directionality. Using the drop down menu above the arrows, the increment of movement can be changed (**Figure 9**), from 0.1 to 1 mm. If necessary, the advanced tab on the left of the page can be opened to provide more information or restore default locations. This tab is generally not needed for autosampler calibration and shouldn't be necessary in most cases.





Figure 9. Set increments.



KinExA Trade-Up Program

Does your lab have an obsolete biosensor that no longer meets your needs? Take advantage of our KinExA Trade-Up Program and upgrade to a KinExA 4000.

Trade-In Value		KinExA 3200: \$ 30,000
		KinExA 3000: \$ 25,000
		SPR or Octet: \$ 20,000
		Autosampler: \$ 10,000

Contract KinExA Research

www.sapidyne.com/cro

With over 25 years of experience you can trust Sapidyne scientists to accurately characterize a wide variety of binding interactions such as antibody-antigen, protein-protein, antibody-cell surface receptors, and more. Our research services are ideal for organizations that require third-party verification, measurement for a special project, or have high affinity interactions and need the most accurate affinity and kinetics data available in the industry. Our research services are also available through Science Exchange or Scientist.com. Standard measurement services include full characterization of a given binding pair using GLP/GMP guidelines, with a detailed formal report provided.

USA

Sapidyne
INSTRUMENTS

info@sapidyne.com

Japan

Scrum Inc
sales@scrum-net.co.jp

Singapore

Yiches pte ltd
sales@yiches.com.sg

Europe

Sapidyne
INSTRUMENTS GMBH

info@sapidyne.eu

Korea

GREENics
greenics@naver.co.kr

China

Beijing Loham Trading Co., Ltd.
Advanced Biotechnologies Limited
sales@advanced-biotechs.com

